

## Oil Free Vacuum Pump



TS15



TS25S



TS45



TS75S



TS75P



TS45S



TS50-V



TS50-P

	TS50-V	TS50-P
<b>MOC</b>	Aluminum Die Cast	Aluminum Die Cast
<b>Vacuum</b>	27"	28"
<b>Flow</b>	50 LPM	50 LPM
<b>Stage</b>	Double Stage	Double Stage
<b>Pressure</b>	110 Psig	110 Psig
<b>Connection</b>	8 mm Hose	8 mm Hose
<b>Weight</b>	10 Kg	10 Kg
<b>Motor</b>	1/4 hp	1/4 hp

	TS15	TS25S	TS45	TS75S	TS75P	TS45S
<b>MOC</b>	Aluminum Die Cast	Aluminum Die Cast	Aluminum Die Cast	Aluminum Die Cast	Aluminum Die Cast	Aluminum Die Cast
<b>Vacuum</b>	22"	27"	22"	27"	22"Hg- 705mm Hg	23"Hg- 550mm Hg
<b>Flow</b>	15 LPM	25 LPM	45 LPM	45 LPM	75 LPM	45 LPM
<b>Stage</b>	Single Stage	Double Stage	Single Stage	Double Stage	Double Stage	Single Stage
<b>Pressure</b>	25 Psig	35 Psig	25 Psig	55 Psig	40 Psig	35 Psig
<b>Connection</b>	8 mm Hose	8 mm Hose	10 mm Hose	10 mm Hose	10 mm Hose	12 mm Hose
<b>Weight</b>	3 Kg	6 Kg	6.5 kg	10 Kg	10 Kg	6.5Kg
<b>Motor</b>	1/20 hp	1/16 hp	1/8 hp	1/4 hp	1/4 hp	1/8 hp
<b>Price</b>						

## Chemical resistant pumps



TS45 CRP



TS75 CRP

	TS45 CRP	TS75 CRP
<b>MOC</b>	Aluminum Die Cast	Aluminum Die Cast
<b>Vacuum</b>	22" Hg	28" Hg
<b>Flow</b>	30 LPM	50 LPM
<b>Stage</b>	Single Stage	Double Stage
<b>Pressure</b>	35 Psig	45 Psig
<b>Connection</b>	12 mm Hose	12 mm Hose
<b>Weight</b>	8Kg	10.5Kg
<b>Motor</b>	1/8 hp	1/4 hp